



Why -8:00 to +6:00 on availability chart

Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	65	6.75	6.77	-5.50	12.2	10.5
+3.25 to +4.50	65	5.75	5.84	-5.50	6.7	5.9
+1.75 to +3.00	70	5.00	4.93	-5.50	6.1	6.5
-2.00 to +1.50	70	4.00	3.99	-7.00	4.6	8.1
-6.00 to -2.25	70	3.00	2.98	-8.75	6.4	12.9
-10.00 to -6.25	70	2.00	2.01	-8.75	10.2	18.0

Prescription Selection Chart



Fitting Height = 14mm min.