


Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+7.25 to +8.00	70	8.50	8.46	-7.75	10.3	8.5
+6.25 to +7.00	70	7.25	7.26	-7.75	9.5	9.6
+4.75 to +6.00	75	6.25	6.25	-5.50	8.3	7.4
+3.25 to +4.50	75	5.50	5.44	-5.50	7.1	7.3
+1.75 to +3.00	75	4.75	4.82	-5.00	8.4	8.0
-2.00 to +1.50	75	4.00	3.97	-5.50	7.0	8.6
-6.00 to -2.25	75	3.00	3.10	-5.50	8.3	11.3
-10.00 to -6.25	75	1.50	1.49	-5.50	8.2	13.0



Fitting Height = 18mm min.

Prescription Selection Chart

