


| Range (D) | Effective Diameter (mm) | Nominal Base | TRUE Curve N=1.530 | Back Curve | Center Thickness (mm) | Edge Thickness (mm) |
|-----------------|-------------------------|--------------|--------------------|------------|-----------------------|---------------------|
| +7.25 to +8.00 | 70 | 8.50 | 8.46 | -7.75 | 10.3 | 8.5 |
| +6.25 to +7.00 | 70 | 7.25 | 7.26 | -7.75 | 9.5 | 9.6 |
| +4.75 to +6.00 | 75 | 6.25 | 6.25 | -5.50 | 8.3 | 7.4 |
| +3.25 to +4.50 | 75 | 5.50 | 5.44 | -5.50 | 7.1 | 7.3 |
| +1.75 to +3.00 | 80 | 4.75 | 4.82 | -5.00 | 8.4 | 8.0 |
| -2.00 to +1.50 | 80 | 4.00 | 3.97 | -5.50 | 7.0 | 8.6 |
| -6.00 to -2.25 | 80 | 3.00 | 3.10 | -5.50 | 8.3 | 11.3 |
| -13.00 to -6.25 | 75 | 1.50 | 1.49 | -5.50 | 8.2 | 13.0 |



Fitting Height = 18mm min.

Prescription Selection Chart

