



Range (D)	Effective Diameter (mm)	Nominal Base	TRUE Curve N=1.530	Back Curve	Center Thickness (mm)	Edge Thickness (mm)
+4.75 to +6.00	75	6.25	6.26	-6.15	10.8	10.7
+3.25 to +4.50	75	5.50	5.50	-6.15	10.6	11.4
+1.75 to +3.00	75	4.75	4.58	-6.15	10.3	12.4
-2.00 to +1.50	75	4.00	4.01	-6.15	9.4	12.3
-6.00 to -2.25	75	3.00	2.88	-4.24	12.3	14.0

Prescription Selection Chart



Fitting Height = 18 mm min.