

Product Reference Guide

1.50 Clear Blue Filter SFSV

Base (mm)	True Curve	Front Curve	Back Radius	Back SAG*	CT (mm)	ET
Nominal	Nominal					
.50	.20	2650.00	-5.60	-92.982	8.5	16.0
1.25	1.00	530.000	-5.60	-92.982	9.8	16.0
2.25	1.98	267.677	-5.60	-92.982	11	16.0
3.25	2.98	177.852	-5.60	-92.982	9.2	13.0
4.25	3.99	132.832	-5.60	-92.982	10.5	13.0
5.25	5.00	106.000	-5.60	-92.982	12.0	13.0
6.25	6.06	88.04	-5.60	-92.982	9.0	8.0
7.25	6.92	76.812	-5.60	-92.982	9.0	7.5
8.25	8.09	65.839	-5.60	-92.982	11	7.0
9.25	8.92	59.417	-5.60	-92.982	13	7.5

*Utilized for power calculation only, not for thickness calculation

