

Product Reference Guide

Plastic

70mm Finished Single Vision Uncoated

Sphere Power Up to -2.00 cyl	Base Curve	Curve Radius	Center Thickness
+4.00	7.00	75.714	6.4
+3.75	6.75	78.519	6.1
+3.50	6.50	81.538	5.6
+3.25	6.25	84.800	5.4
+3.00	7.00	75.714	5.2
+2.75	6.75	78.519	4.9
+2.50	6.50	81.538	4.5
+2.25	6.25	84.800	4.2
+2.00	6.50	81.538	4.0
+1.75	6.25	84.800	3.6
+1.50	6.00	88.333	3.3
+1.25	5.75	92.174	2.9
+1.00	6.00	88.333	2.7
+0.75	5.75	92.174	2.5
+0.50	5.50	96.364	2.5
+0.25	5.25	100.952	2.0
Plano	5.00	106.000	2.0
-0.25	4.75	111.579	2.0
-0.50	4.50	117.778	2.0
-0.75	4.25	124.706	1.9
-1.00	4.00	132.500	1.9
-1.25	3.75	141.333	1.8
-1.50	3.50	151.429	1.8
-1.75	3.25	163.077	1.8
-2.00	3.00	176.667	1.8
-2.25	2.75	192.727	1.7
-2.50	2.50	212.000	1.7
-2.75	2.25	235.556	1.7
-3.00	2.00	265.000	1.7
-3.25	1.75	302.857	1.7
-3.50	1.50	353.333	1.7
-3.75	2.75	192.727	1.7
-4.00	2.50	212.000	1.6
-4.25	2.25	235.556	1.6
-4.50	2.00	265.000	1.6
-4.75	2.25	235.556	1.6
-5.00	2.00	265.000	1.6
-5.25	1.75	302.857	1.6
-5.50	1.50	353.333	1.6
-5.75	1.25	424.000	1.6
-6.00	1.00	530.000	1.6

