



Product Reference Guide



Plastic

Transitions® Signature VII Single Vision Direct Surface Semi-Finished Single Vision (Gray, Brown)

Base Curve	True Curve	Front Radius	50mm SAG*	Back Curve	CT (mm) Nominal	ET (mm) Nominal	Blank Size
0.50	0.44	1217.07	0.26	4.00	10.6	15.8	77
2.00	1.99	265.43	1.18	2.03	16.7	16.8	77
3.50	3.48	152.45	2.06	2.03	18.8	16.7	77
4.25	4.25	124.70	2.53	4.22	14.8	14.8	77
5.25	5.20	102.03	3.11	4.22	11.9	10.3	77
6.25	6.19	85.68	3.73	6.33	13.6	13.7	77
7.25	7.22	73.45	4.39	6.33	13.6	12.0	77
8.75	8.76	60.48	5.41	6.33	12.0	8.4	73
10.5	10.35	51.23	6.51	6.33	14.8	9.5	69

*Utilized for power calculation only, not for thickness calculation.