

Product Reference Guide

Glass

Double Segments

71mm Dual D28 (White, PhotoGray Extra®)

Base Curve	True Curve (1.530)	Front Radius	50mm SAG*	Back Curve	Back Radius	CT (mm) Nominal	ET (mm) Nominal
4.25	4.15	127.711	2.471	-6.00	-88.33	7.5	9.9
6.25	6.20	85.484	3.738	-6.00	-88.33	6.6	6.3

Standard Separation: 13mm

28 Seg

Minimum Separation: 13mm

Base Curve	White	PGX
4.25	+2.00 over +2.00	+2.00 over +2.00
	+2.50 over +2.50	+2.50 over +2.50
6.25	+1.75 over +1.75	+1.75 over +1.75
	+2.00 over +2.00	+2.00 over +2.00
	+2.25 over +2.25	+2.25 over +2.25
	+2.50 over +2.50	+2.50 over +2.50

*Utilized for power calculation only, not for thickness calculation

[Click Here](#) for Additional Multifocal Technical Information

Product Reference Guide

Glass

Double Segments

75mm Dual D35 (PhotoGray Extra®)

Base Curve	True Curve (1.530)	Front Radius	50mm SAG*	Back Curve	Back Radius	CT (mm) Nominal	ET (mm) Nominal
6.25	6.20	85.484	3.738	-6.00	-88.33	7.0	6.7

Standard Separation: 15mm

35 Seg

Minimum Separation: 15mm

Base Curve	PGX
6.25	+2.00 over +2.00 +2.25 over +2.25

*Utilized for power calculation only, not for thickness calculation

[Click Here](#) for Additional Multifocal Technical Information